



Instrument reprocessing
optimally documented
EXPECT MORE

Maintain the overview of all reprocessing procedures and all medical devices in your practice. Prepare checklists and display them at the relevant time for completion. Manage your instruments and filters with an intuitively operated user interface and let us manage tasks and routine tests. And intelligent monitoring in the background eliminates the possibility of operator errors.

Sego® Simple operation, customisable application, secure results.

Function	Sego®4Star***	Sego®4Basic
Fully automatic process documentation	✓	✓
Platform-independent	✓	✓
Multiple workstations	✓	✓
Print labels	✓	✓
VDDS-media interface	✓	✓
User administration with task allocation	✓	
Medical device management	✓	
Device log	✓	
Task management & administration	✓	
Complete instrument management	✓	
Instructions	✓	
Routine tests	✓	
Checklists	✓	
Intelligent SegoAgent	✓	✓
Sego ASSIGN plus* batch assignment only	✓*	✓*
Sego ASSIGN plus* for batch and sterilegoods assignment	✓*	

* for an additional charge



Update, Support & Service
Stay selective

Stay up to date. Updates are simple and quick for you. In the cloud and not in the practice. Our flat-rate service charges are kept simple and leave room for you to develop the practice and to make changes in your devices.

Simply select in **Sego4Basic** or **Sego4Star** the following option independently of the desired functional scope:

- ✓ **mySego.pass** to link up to 4 devices or the extended level
- ✓ **mySego.pass extended** for up to 8 devices
- ✓ **mySego.pass ASSIGN plus** for batch and sterilegoods assignment.

Platform-independent



All brands, logos and trademarks used are the property of the applicable owner or manufacturer.

Sold exclusively by certified dealers.

com|GO|tec®

Gutenbergstr. 3 | D-85716 Unterschleißheim | Phone: +49 89 3270 889-0

mySego®



Instrument reprocessing **optimally** documented

MORE solutions.
MORE quality.
MORE hygiene.



www.sego-soft.info



Medical Device Management & Task Administration

Sego® reminds you

With the wide range of options offered by the mySego technology platform Sego4Star is the most comprehensive and intuitive tool for management of all the medical devices in the practice. All medical devices and instruments can be completely documented and traced. Inventory lists can be generated with one click. The close linkage of task and instrument management enables the display of service intervals and schedules for validation or STK as defined by the user and for



processing. All relevant documents can also be found in the comprehensive archive and the special medical device list.



Routine testing & batch control

Important things at the right time

Sego4Star takes charge of the monitoring of all routine tests and required batch controls. Create additional checklists for all of your devices. MySego displays all tasks that must be completed to the responsible person at the right time.

- ✓ batch-based testing, e.g. Helix
- ✓ cycle-based check, e.g. filter every 20 cycles
- ✓ scheduled check, e.g. seal check, weekly
- ✓ time-critical inspection, e.g. service due




Instrument management & Filter list management

Maintain your overview

You can also manage your instruments and sets with the mySego technology platform, which is closely linked to the task and device management. You can keep a close watch on specialised instruments or instruments and sets that can only be reprocessed to a limited extent and identify ageing inventory. Every Sego4Star user has access to the wide range of functions.

mySego®



Security Relax

Sego offers reliable documentation of the disinfection and sterilization process which can be traced precisely at any time. Recording in combination with electronically signed documents is indispensable for quality assurance and protects you against liability risks. Long-term preservation is ensured by the PDF/A-1 format, which is defined by the ISO 19005-1:2005 standard. The format is electronically signed and as such is tamper-proof and thus legally valid.


Rely on maximum data security. Your data will be in German computer centres that are certified by international standards such as ISO 9001:2008 for quality management systems and ISO/IEC 27001:2013 for systems used for information security management. Compliant with the guideline pursuant to §75b SGB V on the requirements for ensuring IT security.



Support

Sego® Agent

The SegoAgent monitors all activities in the background and only becomes active if discrepancies are detected. It displays dialogues and warning messages based on the situation or offers help with, for example, release decisions.



It displays dialogues and warning messages based on the situation or offers help with, for example, release decisions.